

Specification Sheet



Product Name: DL-AP3

Item No.: MED11722

Purity Specification: 98%

Molecular Formula : C₃H₈NO₅P

CAS Number: 20263-06-3

Formula Weight: 169.073


Storage: Store at RT

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816