

Specification Sheet



Product Name: Okadaic acid(solution)

Item No.: MED11727

Purity Specification: 98%

Molecular Formula : C₄₄H₆₈O₁₃

CASNumber: 78111-17-8

FormulaWeight: 805.003


Storage: Desiccate at -20 °C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816