

Specification Sheet



Product Name: Tasquinimod

Item No.: MED11952

Purity Specification: 98%

Molecular Formula : C₂₀H₁₇F₃N₂O₄

CAS Number: 254964-60-8

Formula Weight: 406.355


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816