

Specification Sheet



Product Name: Losartan Potassium (DuP 753)

Item No.: MED12137

Purity Specification: 98%

Molecular Formula : C₂₂H₂₃CLKN₆O

CASNumber: 124750-99-8

FormulaWeight: 462.009


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816