

# Specification Sheet



**Product Name:** Captopril

*Item No.:* MED12159

*Purity Specification:*  $\geq 98\%$

*Molecular Formula :* C<sub>9</sub>H<sub>15</sub>NO<sub>3</sub>S

*CASNumber:* 62571-86-2

*FormulaWeight:* 217.285


*Storage:* -20 ° C

## Test

## Specification

|        | Min. | Max. |
|--------|------|------|
| PURITY | 98%  | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: [www.med-bio.cn](http://www.med-bio.cn) Tel: 400-186-1816