

Specification Sheet



Product Name: 8-Prenylnaringenin

Item No.: MED12290

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₀H₂₀O₅

CASNumber: 53846-50-7

FormulaWeight: 340.370


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816