

Specification Sheet



Product Name: PACAP 1-27

Item No.: MED12708

Purity Specification: 98%


Molecular Formula : C₁₄₂H₂₂₄N₄₀O₃₉S

CAS Number: 127317-03-7 **Formula Weight:** 3147.608

Storage: Desiccate at -20 °C

| Test | Specification | |
|--------|---------------|------|
| | Min. | Max. |
| PURITY | 98% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816