

Specification Sheet



Product Name: Xenin 8

Item No.: MED12777

Purity Specification: 98%


Molecular Formula : C₅₁H₇₉N₁₅O₉

CASNumber: 117442-28-1 **FormulaWeight:**

Storage: Desiccate at -20 °C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816