

Specification Sheet



Product Name: PG 931

Item No.: MED12937

Purity Specification: $\geq 98\%$

Molecular Formula :

CASNumber: 667430-81-1 FormulaWeight:

Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816