

Specification Sheet



Product Name: Leukotriene B₄

Item No.: MED13137

Purity Specification: ≥95%

Molecular Formula : C₂₀H₃₂O₄

CASNumber: 71160-24-2

FormulaWeight: 336.466


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	95%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816