

Specification Sheet



Product Name: Epoprostenol

Item No.: MED13209

Purity Specification: 98%

Molecular Formula : C₂₀H₃₁O₅

CAS Number: 61849-14-7

Formula Weight: 351.457


Storage: Desiccate at -80 °C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816