

Specification Sheet



Product Name: TC-O 9311

Item No.: MED13300

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₀H₁₉N₃O₄

CASNumber: 444932-31-4 **FormulaWeight:** 365.380

Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816