

Specification Sheet



Product Name: XAP044

Item No.: MED13358

Purity Specification: 98%


Molecular Formula : C₁₅H₉O₄

CASNumber: 196928-50-4 **FormulaWeight:** 380.134

Storage: -20 ° C

| Test | Specification | |
|--------|---------------|------|
| | Min. | Max. |
| PURITY | 98% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816