

Specification Sheet



Product Name: AG-041R

Item No.: MED13397

Purity Specification: $\geq 98\%$


Molecular Formula : C₃₁H₃₆N₄O₅

CASNumber: 159883-95-1 **FormulaWeight:** 544.641

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816