

Specification Sheet



Product Name: O-Phospho-D-Serine

Item No.: MED13434

Purity Specification: $\geq 98\%$

Molecular Formula : C₃H₈NO₆P

CAS Number: 73913-63-0

Formula Weight: 185.073


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816