

Specification Sheet



Product Name: 5,7-Dichlorokynurenic acid

Item No.: MED14317

Purity Specification: 98%


Molecular Formula : C₁₀H₅CL₂NO₃

CASNumber: 131123-76-7 **FormulaWeight:** 258.058

Storage: Store at RT

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816