

Specification Sheet



Product Name: Paxilline

Item No.: MED14470

Purity Specification: 98%

Molecular Formula : C₂₇H₃₃NO₄

CASNumber: 57186-25-1

FormulaWeight: 435.555


Storage: Desiccate at -20 °C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816