

Specification Sheet



Product Name: Tipifarnib (Zarnestra)

Item No.: MED14745

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₇H₂₂Cl₂N₄O

CASNumber: 192185-72-1 **FormulaWeight:** 489.396


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816