

Specification Sheet



Product Name: Tenofovir Disoproxil Fumarate

Item No.: MED15145

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₃H₃₄N₅O₁₄P

CASNumber: 202138-50-9

FormulaWeight: 635.515


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816