

Specification Sheet



Product Name: Psoralen

Item No.: MED15747

Purity Specification: $\geq 98\%$

Molecular Formula : C₁₁H₆O₃

CASNumber: 66-97-7

FormulaWeight: 186.164


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816