

Specification Sheet



Product Name: Ginsenoside Rg3

Item No.: MED15847

Purity Specification: $\geq 98\%$

Molecular Formula : C₄₂H₇₂O₁₃

CASNumber: 14197-60-5

FormulaWeight: 785.013


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816