

Specification Sheet



Product Name: Isochlorogenic acid A

Item No.: MED15924

Purity Specification: 98%

Molecular Formula : C₂₅H₂₄O₁₂

CAS Number: 2450-53-5

Formula Weight: 516.451


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816