

Specification Sheet



Product Name: Pseudoginsenoside- F11

Item No.: MED16013

Purity Specification: $\geq 98\%$

Molecular Formula : C₄₁H₇₀O₁₄

CASNumber: 69884-00-0

FormulaWeight: 786.986


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816