

# Specification Sheet



**Product Name:** Senkyunolide H

*Item No.:* MED16380

**Purity Specification:** 98%

**Molecular Formula :** C<sub>15</sub>H<sub>11</sub>NO<sub>5</sub>

**CASNumber:** 94596-27-7

**FormulaWeight:** 285.252


**Storage:** -20 ° C

## Test

## Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816