

# Specification Sheet



**Product Name:** Tropisetron

*Item No.:* MED16784

**Purity Specification:** 98%

**Molecular Formula :** C<sub>17</sub>H<sub>21</sub>CLN<sub>2</sub>O<sub>2</sub>

**CASNumber:** 105826-92-4

**FormulaWeight:** 320.814


**Storage:** -20 ° C

**Test**

**Specification**

	<b>Min.</b>	<b>Max.</b>
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816