

Specification Sheet



Product Name: PEAQX

Item No.: MED16810

Purity Specification: 98%


Molecular Formula : C₁₇H₁₇BRN₃O₅P

CASNumber: 459836-30-7 **FormulaWeight:** 454.212

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816