

Specification Sheet



Product Name: O-Phospho-L-serine

Item No.: MED16907

Purity Specification: 98%

Molecular Formula : C₃H₆NO₆P

CASNumber: 407-41-0

FormulaWeight: 183.057


Storage: Store at RT

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816