

Specification Sheet



Product Name: trans-Triprolidine hydrochloride

Item No.: MED16950

Purity Specification: 98%

Molecular Formula : C₁₉H₂₃CLN₂


CASNumber: 550-70-9

FormulaWeight: 314.852

Storage: Store at RT

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816