

Specification Sheet



Product Name: Tonabersat

Item No.: MED17368

Purity Specification: 98%


Molecular Formula : C₂₀H₁₉ClFNO₄

CASNumber: 175013-84-0 **FormulaWeight:** 391.821

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816