

Specification Sheet



Product Name: Betamethasone Dipropionate

Item No.: MED17398

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₇H₃₄ClFO₇

CASNumber: 5593-20-4

FormulaWeight: 525.006


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816