

# Specification Sheet



**Product Name:** Dipiperodon HCl

*Item No.:* MED17409

**Purity Specification:** 98%

**Molecular Formula :** C<sub>22</sub>H<sub>28</sub>CLN<sub>3</sub>O<sub>4</sub>


**CASNumber:** 537-12-2

**FormulaWeight:** 433.928

**Storage:** -20 ° C

| Test   | Specification |      |
|--------|---------------|------|
|        | Min.          | Max. |
| PURITY | 98%           | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816