

Specification Sheet



Product Name: L-Kynurenine

Item No.: MED17453

Purity Specification: $\geq 98\%$

Molecular Formula : C₁₀H₁₂N₂O₃

CASNumber: 2922-83-0

FormulaWeight: 208.214


Storage: Store at RT

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816