

Specification Sheet



Product Name: LXR-623

Item No.: MED17483

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₁H₁₂ClF₅N₂

CASNumber: 875787-07-8 **FormulaWeight:** 422.778


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816