

Specification Sheet



Product Name: Palosuran

Item No.: MED17493


Purity Specification: $\geq 98\%$

Molecular Formula : C₂₅H₃₀N₄O₂

CASNumber: 540769-28-6 FormulaWeight: 418.531 Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816