

Specification Sheet



Product Name: XL335

Item No.: MED17503

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₅H₂₄F₂N₂O₃

CAS Number: 629664-81-9 **Formula Weight:** 438.466


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816