

Specification Sheet



Product Name: Aliskiren Hemifumarate

Item No.: MED17532

Purity Specification: $\geq 98\%$

Molecular Formula : C₆₄H₁₁₀N₆O₁₆

CAS Number: 173334-58-2 **Formula Weight:** 1219.589


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816