

Specification Sheet



Product Name: AGN 194310

Item No.: MED17574

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₈H₂₄O₂S

CASNumber: 229961-45-9 **FormulaWeight:** 424.554


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816