

Specification Sheet



Product Name: MNI caged kainic acid

Item No.: MED17773

Purity Specification: $\geq 98\%$

Molecular Formula :

CASNumber: 1315378-75-6 *FormulaWeight:*


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816