

Specification Sheet



Product Name: Ethacrynic acid - d5

Item No.: MED17853

Purity Specification: $\geq 99\%$


Molecular Formula : C₁₃H₁₂CL₂O₄

CASNumber: 1330052-59-9 **FormulaWeight:** 308.169

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	99%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816