

Specification Sheet



Product Name: SynaptoRed™ C2

Item No.: MED17859

Purity Specification: 98%

Molecular Formula : C₃₀H₄₅BR₂N₃

CAS Number: 162112-35-8

Formula Weight: 607.506


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816