

Specification Sheet



Product Name: Iberin(solution)

Item No.: MED17891

Purity Specification: $\geq 95\%$

Molecular Formula : C₅H₉NOS₂

CASNumber: 505-44-2

FormulaWeight: 163.261


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	95%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816