

Specification Sheet



Product Name: 4-Dafachronic Acid

Item No.: MED17942

Purity Specification: 98%

Molecular Formula : C₂₇H₄₂O₃

CASNumber: 23017-97-2

FormulaWeight: 414.621


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816