

Specification Sheet



Product Name: PI4KIII beta inhibitor 3

Item No.: MED18111

Purity Specification: 98%


Molecular Formula : C₂₂H₂₂N₈O₅

CASNumber: 1245319-54-3 **FormulaWeight:** 446.528

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816