

Specification Sheet



Product Name: BX517(PDK1 inhibitor2)

Item No.: MED18182

Purity Specification: 98%

Molecular Formula : C₁₅H₁₄N₄O₂

CASNumber: 850717-64-5 **FormulaWeight:** 282.297

Storage: -20 ° C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816