

Specification Sheet



Product Name: Aprotinin

Item No.: MED18397

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₈₄H₄₃₂N₈₄O₇₉S₇

CAS Number: 9087-70-1

Formula Weight: 6511.439


Storage: Store at 2-8 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816