

Specification Sheet



Product Name: Cathepsin S inhibitor

Item No.: MED18413

Purity Specification: 98%


Molecular Formula : C₂₄H₃₀N₂O₆

CASNumber: 1373215-15-6 **FormulaWeight:** 442.505

Storage: -20 ° C

| Test | Specification | |
|--------|---------------|------|
| | Min. | Max. |
| PURITY | 98% | 100% |

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816