

Specification Sheet



Product Name: Saquinavir mesylate

Item No.: MED18432


Purity Specification: $\geq 98\%$

Molecular Formula : C₃₉H₅₄N₆O₈S

CASNumber: 149845-06-7 FormulaWeight: 766.946 Storage: -20 °C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816