

Specification Sheet



Product Name: Phosphoramidon Disodium Salt

Item No.: MED18519

Purity Specification: 98%


Molecular Formula : C₂₃H₃₄N₃Na₂O₁₀P

CAS Number: 164204-38-0 **Formula Weight:** 589.483

Storage: Desiccate at -20 °C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816