

Specification Sheet



Product Name: Hexa-D-arginine

Item No.: MED18529

Purity Specification: 98%

Molecular Formula : C₃₆H₇₅N₂₅O₆

CASNumber: 673202-67-0

FormulaWeight: 954.145


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816