

Specification Sheet



Product Name: KNK437

Item No.: MED18531

Purity Specification: $\geq 98\%$

Molecular Formula : C₁₃H₁₁NO₄

CASNumber: 218924-25-5 **FormulaWeight:** 245.231


Storage: -20 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816